

# 96-Well Plasmid Kit

*For research use only*

<b>Sample</b>	: 1-2 ml of cultured bacterial cells
<b>Binding Capacity</b>	: 10 µg/well
<b>Operation</b>	: centrifuge/vacuum manifold
<b>Yield</b>	: 5-10 µg for high-copy plasmid 0.5-5 µg for low-copy plasmid
<b>Operation time</b>	: within 60 minutes
<b>Application</b>	: Fluorescent/Radioactive Sequencing, Restriction Enzyme Digestion, Library Screening, Ligation, Transformation and PCR

**Geneaid**



CERTIFICATE NO. 6501225407  
[www.geneaid.com](http://www.geneaid.com)

## Introduction

The 96-Well Plasmid Kit was designed for rapid isolation of plasmid DNA from 1-2 ml of cultured bacterial cells. Modified Alkaline Lysis method (1) and RNase A treatment are used to obtain clear cell lysate with minimal genomic DNA/RNA contaminants. In the presence of chaotropic salt, plasmid DNA in the lysate binds to the glass fiber matrix (2) of each well of the plate. Contaminants are removed with a Wash Buffer (containing ethanol) and the purified plasmid DNA is eluted by a low salt Elution Buffer or TE. Typical yields are 5-10 µg for high-copy number plasmid or 0.5-5 µg for low-copy number plasmid. DNA phenol extraction or alcohol precipitation is not required and the entire procedure can be completed within 1 hour. The purified plasmid DNA is ready for use in various downstream applications.

## Quality Control

The quality of the 96-Well Plasmid Kit is tested on a lot-to-lot basis by isolating plasmid DNA from a 1.5 ml overnight *E. coli* (DH5α) culture, containing plasmid pBluescript (A600 > 2 U/ml). More than 8 µg of plasmid DNA can be quantified with a spectrophotometer. The purified plasmid (1 µg) is used in EcoR I digestion, and checked by electrophoresis.

## Kit Contents

Name	PDA02	PDA04	PDA10
PD 1 Buffer*	25 ml	65 ml	65 ml x 2
PD 2 Buffer**	25 ml	75 ml	75 ml x 2
PD 3 Buffer	45 ml	100 ml	100 ml x 2
W1 Buffer	60 ml	130 ml	130 ml x 2
Wash Buffer*** (Add Ethanol)	12.5 ml (50 ml)	25 ml (100 ml)	50 ml (200 ml)
Elution Buffer (10 mM Tris-HCl, pH 8.5 at 25°C)	30 ml	60 ml	120 ml
RNase A (50 mg/ml)	50 µl	130 µl	130 µl x 2
Plasmid Plate	2 pcs	4 pcs	10 pcs
0.35 ml Collection Plate	2 pcs	4 pcs	10 pcs
Adhesive Film	4 pcs	8 pcs	20 pcs

## Order Information

Product Name	Package Size	Cat. No.
High-Speed Plasmid Mini Kit	100/300 preps	PD100/300
High-Speed Plasmid Mini Kit (10-50 Kb)	100/300 preps	PDL100/300
High-Speed Plasmid Advance Kit (50-100 ml bacterial culture)	25 preps	PA025
Geneaid Plasmid Midi Kit	25 preps	PI025
Geneaid Plasmid Midi Kit (Endotoxin Free)	25 preps	PIE25
Geneaid Plasmid Maxi Kit	10/25 preps	PM010/025
Geneaid Plasmid Maxi Kit (Endotoxin Free)	10/25 preps	PME10/E25
96-Well Plasmid Kit	2/4/10 x 96 Wells	PDA02/04/10
Vacuum Manifold (Accessories)	1 Set	ZVF01

\*Add provided RNase A to the PD1 Buffer and store at 4°C

\*\*If precipitates have formed in the PD2 Buffer, warm the buffer in a 37°C water bath to dissolve

\*\*\*Add absolute ethanol (see the bottle label for volume) to the Wash Buffer prior to initial use

## Caution

PD3 Buffer and W1 Buffer contain guanidine hydrochloride which is a harmful irritant. During the procedure, always wear a lab coat, disposable gloves, and protective goggles.

## References

- (1) Birnboim, H. C., and Doly, J. (1979) *Nucleic Acids Res.* 7, 1513.
- (2) Vogelstein, B. and Gillespie, D. (1979) *Proc. Natl. Acad. Sci. USA* 76, 615

## 96-Well Plasmid Kit Centrifuge Protocol

- Additional requirements: centrifugation system for 96-well plates, 2 ml collection plates, absolute ethanol

Step 1 Harvesting	<ul style="list-style-type: none"> <li>Transfer 1.5 ml of cultured bacterial cells to a 2 ml collection plate and centrifuge for 10 minutes at 3,000 x g.</li> <li>Remove the culture medium.</li> </ul>
Step 2 Re-suspension	<ul style="list-style-type: none"> <li>Add <b>100 µl of PD1 Buffer</b> (RNase A added) to each well of the plate and re-suspend the cell pellet by pipetting or vortex.</li> </ul>
Step 3 Lysis	<ul style="list-style-type: none"> <li>Add <b>100 µl of PD2 Buffer</b> to each well and mix gently by shaking.</li> <li>Let stand at room temperature for at least 2 minutes to ensure the lysate is clear.</li> </ul>
Step 4 Neutralization	<ul style="list-style-type: none"> <li>Add <b>150 µl of PD3 Buffer</b> and mix immediately by shaking gently, then centrifuge for 10 minutes at 3,000 x g. Following centrifugation, a denatured protein pellet and clear lysate (supernatant) will be present in each well.</li> </ul>
Step 5 DNA Binding	<ul style="list-style-type: none"> <li>Place a <b>Plasmid Plate</b> on a new 2 ml collection plate.</li> <li>Carefully transfer <b>ONLY</b> the clear lysate (supernatant) from Step 4 to the <b>Plasmid Plate</b>.</li> <li>Centrifuge for 5 minutes at 3,000 x g.</li> <li>Discard the flow-through and place the <b>Plasmid Plate</b> back on the 2 ml collection plate.</li> </ul>
Step 6 Wash	<ul style="list-style-type: none"> <li>Add <b>250 µl of W1 Buffer</b> into each well of the <b>Plasmid Plate</b> and centrifuge for 5 minutes at 3,000 x g.</li> <li>Add <b>250 µl of Wash Buffer</b> (ethanol added) into each well of the <b>Plasmid Plate</b> to wash again and centrifuge for 5 minutes at 3,000 x g.</li> <li>Discard the flow-through and place the <b>Plasmid Plate</b> back on the 2 ml collection plate.</li> <li>Centrifuge for 10 minutes at 3,000 x g to remove any ethanol residue.</li> </ul>
Step 7 DNA Elution	<ul style="list-style-type: none"> <li>Transfer the <b>Plasmid Plate</b> to a <b>0.35 ml Collection Plate</b>.</li> <li>Add <b>50-100 µl of Elution Buffer</b> or TE into the center of each membrane matrix.</li> <li>Let stand for at least 2 minutes to ensure the <b>Elution Buffer</b> or TE is absorbed.</li> <li>Centrifuge for 5 minutes at 3,000 x g to elute the purified DNA.</li> </ul>

## 96-Well Plasmid Kit Vacuum Protocol

- Additional requirements: centrifuge system for 96-well plates, multi-well plate vacuum manifold, 2 ml collection plate, absolute ethanol

Step 1 Harvesting	<ul style="list-style-type: none"> <li>Transfer 1.5 ml of cultured bacterial cells to a 2 ml collection plate and centrifuge for 10 minutes at 3,000 x g.</li> <li>Remove the supernatant.</li> </ul>
Step 2 Re-suspension	<ul style="list-style-type: none"> <li>Add <b>100 µl of PD1 Buffer</b> (RNase A added) to each well of the plate and re-suspend the cell pellet by pipetting or vortex.</li> </ul>
Step 3 Lysis	<ul style="list-style-type: none"> <li>Add <b>100 µl of PD2 Buffer</b> and mix gently by shaking.</li> <li>Let stand at room temperature for at least 2 minutes to ensure the lysate is clear.</li> </ul>
Step 4 Neutralization	<ul style="list-style-type: none"> <li>Add <b>150 µl of PD3 Buffer</b> and mix immediately by shaking gently and centrifuge for 10 minutes at 3,000 x g. Following centrifugation, a denatured protein pellet and clear lysate (supernatant) will be present in each well.</li> </ul>
Step 5 DNA Binding	<ul style="list-style-type: none"> <li>Place a 2 ml collection plate on the base of the vacuum manifold and place a <b>Plasmid Plate</b> on top.</li> <li>Carefully transfer <b>ONLY</b> the clear lysate (supernatant) from Step 4 to the <b>Plasmid Plate</b> (approximately 350 µl).</li> <li>Apply vacuum at 10 inches Hg for at least 5 minutes to ensure the wells have emptied.</li> </ul>
Step 6 Wash	<ul style="list-style-type: none"> <li>Turn off the vacuum pump and add <b>250 µl of W1 Buffer</b> to each well of the <b>Plasmid Plate</b>.</li> <li>Apply vacuum at 10 inches Hg for at least 5 minutes to ensure the wells have emptied.</li> <li>Turn off the vacuum pump and add <b>250 µl of Wash Buffer</b> (ethanol added) to each well of the <b>Plasmid Plate</b>.</li> <li>Apply vacuum at 10 inches Hg for at least 5 minutes to ensure the wells have emptied.</li> <li>Apply vacuum for an additional 10 minutes (or incubate at 60°C for 5-10 minutes) to remove any ethanol residue.</li> </ul>
Step 7 DNA Elution	<ul style="list-style-type: none"> <li>Turn off the vacuum pump and transfer the <b>Plasmid Plate</b> to a <b>0.35 ml Collection Plate</b>.</li> <li>Add <b>50-100 µl of Elution Buffer</b> or TE into the center of each membrane matrix.</li> <li>Let stand for at least 2 minutes to ensure the <b>Elution Buffer</b> or TE is absorbed.</li> <li>Centrifuge for 5 minutes at 3,000 x g to elute the purified DNA.</li> </ul>

## Troubleshooting

Problem	Possible Reasons/Solution
Low Yield	<p><b>Bacterial cells were not lysed completely</b></p> <ul style="list-style-type: none"> <li>If more than 10 OD<sub>600</sub> units of bacterial culture are used, dilute into multiple tubes.</li> <li>Following PD3 Buffer addition, pipetting or gentle shaking will help to ensure the sample is homologous.</li> </ul> <p><b>Incorrect DNA Elution Step</b></p> <ul style="list-style-type: none"> <li>Ensure that Elution Buffer is added into the center of each well.</li> </ul> <p><b>Incomplete DNA Elution</b></p> <ul style="list-style-type: none"> <li>If plasmid DNA are larger than 10 Kb, use pre-heated Elution Buffer (60~70°C) in the Elution step.</li> </ul>
Eluted DNA does not perform well in downstream applications	<p><b>Residual Ethanol contamination</b></p> <ul style="list-style-type: none"> <li>Following Wash step, dry the Plasmid Plate with extended vacuum/centrifugation time of approximately 5 minutes.</li> </ul> <p><b>RNA contamination</b></p> <ul style="list-style-type: none"> <li>Prior to using PD1 Buffer, ensure RNase A is added.</li> </ul> <p><b>Genomic DNA contamination</b></p> <ul style="list-style-type: none"> <li>Do not use overgrown bacterial cultures.</li> <li>During PD2 and PD3 Buffer addition, mix gently to prevent genomic DNA shearing.</li> </ul>